

Draining rack

- Drying racks makes it easy to dry glassware in the laboratory.
- Our range includes draining racks with both fixed and changeable rods.
- Model ETTA is made entirely of polypropylene. These racks are also available in tailored sizes. Ask more from our experts!
- The material of model ETTB draining rack is PVC, rods polypropylene. The rods are individually placed into the holes of the draining rack and are clicked into place. Their positions can be changed at any time.




Product information

| Product code and standard dimensions width-height [mm] | Model | Additional information |
|--|------------------------------|--|
| ETTA-500-700 | A, fixed rods | Containment sump available on request. |
| ETTA-600-700 | A, fixed rods | Containment sump available on request. |
| ETTA-1000-700 | A, fixed rods | Containment sump available on request. |
| ETTB-500-500 | B, changeable rods | Containment sump included. |
| ETTB-600-600 | B, changeable rods | Containment sump included. |
| ETTC-280-312-450 (depth 312) | C, worktop model, fixed rods | Containment sump included. |


Additional information

| Model | Rod lengths [mm] |
|---|--|
| ETTA (A-model, fixed rods) | 190 bottom row, 140 others (as standard) |
| ETTB (B-model, changeable rods) | 100, 60 |
| ETTC (C-model, worktop model, fixed rods) | 95 bottom row, 150 others |

VISU KALUSTE OY

 Eeliksentie 6, 34870 Visuvesi, Finland

 visu@visu.fi

 +358 44 749 4940

 visu.fi